


ClarusScope™



Endoscopic Shunt Placement System



1995
Clarus invents the
NeuroPen




2009
Clarus Invents the
ClarusScope


See the Difference:

- **Reduces costs up to 50%**
- **Easier Set Up**
- **Always in Focus**
- **Superior Optics**
- **No Optical Coupler Needed**

Endoscopic Shunt Placement System



**Reusable ClarusScope
#3000N**



**Disposable ShuntPen
Cannula #2120C
(sold in boxes of 5)**



- **No Optical Coupler is required**
- **ClarusScope is always in focus**
- **Functional dimensions are identical to the NeuroPen**

STERILIZATION: The reusable ClarusScope is compatible with ETO, STERRAD® and Steris® sterilization methods. The single use cannula is supplied sterile.



13355 10th AVE N
Suit 110
Minneapolis, MN 55441

763-525-8403
800-359-2372
FAX: 763-525-8656

Email: admin@clarus-medical.com
Web Site: www.clarus-medical.com